



# In-Situ Characterization of Dry Microstructured Surfaces for Human Skin Adhesion


Aly Hassan<sup>1</sup>, Benjamin Hatton<sup>2</sup>, Tobin Fillete<sup>1</sup>

<sup>1</sup>Department of Mechanical and Industrial Engineering, University of Toronto

<sup>2</sup>Department of Materials Science and Engineering, University of Toronto

## Background

As the use of wearable sensors in a more personalized healthcare system increases, the need for a reliable, reusable, and clean adhesive mechanism to skin is critical.




- Skin adhesive patches enable<sup>[1,2]</sup>:
- Monitoring vital signs (EEG, ECG, EMG)
  - Body motion & vibration sensing
  - Controlled drug delivery

Current state of the art in bioinspired dry adhesive for skin vary from mushroom-shaped pillars to octopus-inspired suction cups. The highest reported normal adhesion to skin is  $\sim 2 \text{ N/cm}^2$ .




Mushroom-shaped pillar array by Kwak et al. in 2011<sup>[3]</sup>  
Adhesion:  $1.3 \text{ N/cm}^2$



Suction cup array adhesive by Baik et al. in 2018<sup>[4]</sup>  
Adhesion:  $1.8 \text{ N/cm}^2$

## Approach




To improve upon dry microstructured adhesion to skin, the mechanisms of adhesion are first understood on commonly studied surfaces (i.e. glass) and then skin-replica surfaces.



The dry microstructured surfaces show a significant drop-off in adhesion when moving from a hard, smooth glass surface to a soft, rough surface (skin replica).

## Results

Commercial mushroom-shaped pillar arrays (nanogriptech, tip diameter  $150\mu\text{m}$ , AR  $\sim 2$ ) were tested on glass surfaces and showed normal adhesion of up to  $13 \text{ N/cm}^2$  with significant dependence on preload. Normal adhesion testing on samples of varying diameters revealed a size-scale effect: smaller samples had higher adhesion per unit area. In-situ characterization was also enabled through frustrated-TIR<sup>[5]</sup>.



Preliminary results reveal a mostly random pillar detachment sequence and a change in the contact behavior at very high preloads.

## Conclusions & Next Steps

- Preliminary results show a random detachment mechanism in most cases, but this needs further statistical analysis.
- Contact appears to change under high preloads, however the nature of this drop in adhesion remains unclear.
- A size-scale effect was observed where smaller samples achieved higher adhesion per unit area.
- Future work will include understanding changes in the contact mechanism with size-scale effect and on skin replica samples.

[1] S. Baik, J. Kim, H.J. Lee, T.H. Lee, C. Pang, Advanced Science 5 (2018) 1870045.

[2] S. Baik, H.J. Lee, D.W. Kim, J.W. Kim, Y. Lee, C. Pang, Advanced Materials (2019) 1803309.

[3] I. Hwang, H.N. Kim, M. Seong, S.-H. Lee, M. Kang, H. Yi, W.G. Bae, M.K. Kwak, H.E. Jeong, Advanced Healthcare Materials 7 (2018) 1800275.

[4] M.K. Kwak, H.-E. Jeong, K.Y. Suh, Advanced Materials 23 (2011) 3949–3953.

[5] V. Tinnemann, L. Hernández, S.C.L. Fischer, E. Arzt, R. Bennewitz, R. Hensel, Advanced Functional Materials 29 (2019) 1807713.